

## Antibodies Product Datasheet

|                               |                                                                                                                                                                             |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CLIC1 Antibody                                                                                                                                                              |
| <b>Catalog Number:</b>        | C47495                                                                                                                                                                      |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                           |
| <b>Applications:</b>          | IHC-p, IF(paraffin section)                                                                                                                                                 |
| <b>Dilutions:</b>             | IHC 1:100-200                                                                                                                                                               |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                            |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.                                                               |
| <b>Immunogen:</b>             | Synthetic Peptide of CLIC1                                                                                                                                                  |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                     |
| <b>Gene Name:</b>             | CLIC1                                                                                                                                                                       |
| <b>Protein Name:</b>          | Chloride intracellular channel protein 1 (Chloride channel ABP) (Nuclear chloride ion channel 27) (NCC27) (Regulatory nuclear chloride ion channel protein) (hRNCC)         |
| <b>Other Name:</b>            | Chloride intracellular channel protein 1 (Chloride channel ABP) (Nuclear chloride ion channel 27) (NCC27) (Regulatory nuclear chloride ion channel protein) (hRNCC); YM3542 |
| <b>Human Swiss Prot No.:</b>  | O00299                                                                                                                                                                      |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                     |
| <b>Specificity:</b>           | CLIC1 protein(A216) detects endogenous levels of CLIC1                                                                                                                      |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                |