

Antibodies Product Datasheet

Product Name:	eEF1A2 binding protein Antibody
Catalog Number:	C47498
Reactivity:	Human, Rat, MousePig
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:500-1:2000 IHC: 1:50-1:200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of eEF1A2 binding protein
Concentration:	1 mg/ml
Gene Name:	IGFN1
Protein Name:	Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1) (KY-interacting protein 1)
Other Name:	Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1) (KY-interacting protein 1); YM3473
Human Swiss Prot No.:	Q86VF2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous eEF1A2 binding protein protein
Storage and Stability:	-20°C/1 year