

Antibodies Product Datasheet

Product Name:	GRIK2 (GluR6) Antibody
Catalog Number:	C47503
Reactivity:	Human
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GRIK2 (GluR6)
Concentration:	1 mg/ml
Gene Name:	GRIK2
Protein Name:	Glutamate receptor, ionotropic kainate 2 (Excitatory amino acid receptor 4) (EAA4) (Glutamate receptor 6) (GluR-6) (GluR6)
Other Name:	Glutamate receptor, ionotropic kainate 2 (Excitatory amino acid receptor 4) (EAA4) (Glutamate receptor 6) (GluR-6) (GluR6); YM3562
Human Swiss Prot No.:	Q13002
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRIK2(GluR6) protein(A240) detects endogenous levels of GRIK2(GluR6)
Storage and Stability:	-20°C/1 year