

Antibodies Product Datasheet

Product Name:	HAO1 Antibody
Catalog Number:	C47504
Reactivity:	Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of HAO1
Concentration:	
Gene Name:	HAO1
Protein Name:	Hydroxyacid oxidase 1
Other Name:	Hydroxyacid oxidase 1; HAOX1; Glycolate oxidase; GOX; YM3437
Human Swiss Prot No.:	Q9UJM8
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous HAO1 protein.
Storage and Stability:	-20°C/1 year