

Antibodies Product Datasheet

Product Name:	IGFN1 Antibody
Catalog Number:	C47514
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:2000-5000 IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of IGFN1
Concentration:	
Gene Name:	IGFN1
Protein Name:	Immunoglobulin-like and fibronectin type III domain-containing protein 1
Other Name:	Immunoglobulin-like and fibronectin type III domain-containing protein 1; EEF1A2-binding protein 1; KY-interacting protein 1; YM3401
Human Swiss Prot No.:	Q86VF2
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous IGFN1 protein.
Storage and Stability:	-20°C/1 year