

Antibodies Product Datasheet

Product Name:	KCNK10 (TREK-2) Antibody
Catalog Number:	C47516
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KCNK10 (TREK-2)
Concentration:	1 mg/ml
Gene Name:	KCNK10
Protein Name:	Potassium channel subfamily K member 10 (Outward rectifying potassium channel protein TREK-2) (TREK-2 K(+) channel subunit)
Other Name:	Potassium channel subfamily K member 10 (Outward rectifying potassium channel protein TREK-2) (TREK-2 K(+) channel subunit); YM3559
Human Swiss Prot No.:	P57789
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNK10(TREK-2) protein(A237) detects endogenous levels of KCNK10(TREK-2)
Storage and Stability:	-20°C/1 year