

Antibodies Product Datasheet

Product Name:	KCNK9 (TASK-3) Antibody
Catalog Number:	C47517
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KCNK9 (TASK-3)
Concentration:	1 mg/ml
Gene Name:	KCNK9
Protein Name:	Potassium channel subfamily K member 9 (Acid-sensitive potassium channel protein TASK-3) (TWIK-related acid-sensitive K(+) channel 3) (Two pore potassium channel KT3.2) (Two pore K(+) channel KT3.2)
Other Name:	Potassium channel subfamily K member 9 (Acid-sensitive potassium channel protein TASK-3) (TWIK-related acid-sensitive K(+) channel 3) (Two pore potassium channel KT3.2) (Two pore K(+) channel KT3.2); YM3561
Human Swiss Prot No.:	Q9NPC2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNK9(TASK-3) protein(A239) detects endogenous levels of KCNK9(TASK-3)
Storage and Stability:	-20°C/1 year