

Antibodies Product Datasheet

Product Name:	KCNN3 (SK3) Antibody
Catalog Number:	C47519
Reactivity:	Human, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KCNN3 (SK3)
Concentration:	1 mg/ml
Gene Name:	KCNN3
Protein Name:	Small conductance calcium-activated potassium channel protein 3 (SK3) (SKCa 3) (SKCa3) (KCa2.3)
Other Name:	Small conductance calcium-activated potassium channel protein 3 (SK3) (SKCa 3) (SKCa3) (KCa2.3); YM3566
Human Swiss Prot No.:	Q9UGI6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNN3(SK3) protein(A246) detects endogenous levels of KCNN3(SK3)
Storage and Stability:	-20°C/1 year