

Antibodies Product Datasheet

Product Name:	Kv10.2 Antibody
Catalog Number:	C47521
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Kv10.2
Concentration:	1 mg/ml
Gene Name:	KCNH5
Protein Name:	Potassium voltage-gated channel subfamily H member 5 (Ether-a-go-go potassium channel 2) (hEAG2) (Voltage-gated potassium channel subunit Kv10.2)
Other Name:	Potassium voltage-gated channel subfamily H member 5 (Ether-a-go-go potassium channel 2) (hEAG2) (Voltage-gated potassium channel subunit Kv10.2); YM3577
Human Swiss Prot No.:	Q8NCM2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Kv10.2 protein(A260) detects endogenous levels of Kv10.2
Storage and Stability:	-20°C/1 year