

Antibodies Product Datasheet

Product Name:	KV β 1 Antibody
Catalog Number:	C47522
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KV β 1
Concentration:	1 mg/ml
Gene Name:	KCNAB1
Protein Name:	Voltage-gated potassium channel subunit beta-1 (K(+)) channel subunit beta-1) (Kv-beta-1)
Other Name:	Voltage-gated potassium channel subunit beta-1 (K(+)) channel subunit beta-1) (Kv-beta-1); YM3572
Human Swiss Prot No.:	Q14722
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KV β 1 protein(A255) detects endogenous levels of KV β 1
Storage and Stability:	-20°C/1 year