

Antibodies Product Datasheet

Product Name:	LDHD Antibody
Catalog Number:	C47523
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of LDHD
Concentration:	
Gene Name:	LDHD
Protein Name:	Probable D-lactate dehydrogenase, mitochondrial
Other Name:	Probable D-lactate dehydrogenase, mitochondrial (DLD) (Lactate dehydrogenase D) (EC 1.1.2.4); YM3324
Human Swiss Prot No.:	Q86WU2
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous LDHD protein.
Storage and Stability:	-20°C/1 year