

Antibodies Product Datasheet

Product Name:	LexA DNA Binding Region Antibody
Catalog Number:	C47524
Reactivity:	
Applications:	WB
Dilutions:	WB: 1:10000-50000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of LexA DNA Binding Region
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	YM3412
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous LexA DNA Binding Region protein.
Storage and Stability:	-20°C/1 year