

Antibodies Product Datasheet

Product Name:	A-FABP Antibody
Catalog Number:	C47558
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human A-FABP
Concentration:	1 mg/ml
Gene Name:	FABP4
Protein Name:	A-FABP
Other Name:	Fatty acid-binding protein, adipocyte (Adipocyte lipid-binding protein) (ALBP) (Adipocyte-type fatty acid-binding protein) (A-FABP) (AFABP) (Fatty acid-binding protein 4)
Human Swiss Prot No.:	P15090
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human A-FABP
Storage and Stability:	-20°C/1 year