

Antibodies Product Datasheet

Product Name:	Amphiphysin II Antibody
Catalog Number:	C47559
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human Amphiphysin II.
Concentration:	1 mg/ml
Gene Name:	BIN1
Protein Name:	Myc box-dependent-interacting protein 1
Other Name:	BIN1; AMPHL; Myc box-dependent-interacting protein 1; Amphiphysin II; Amphiphysin-like protein; Box-dependent myc-interacting protein 1; Bridging integrator 1
Human Swiss Prot No.:	O00499
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Amphiphysin II Polyclonal Antibody detects endogenous levels of Amphiphysin II protein.
Storage and Stability:	-20°C/1 year