

Antibodies Product Datasheet

Product Name:	Ankyrin B Antibody
Catalog Number:	C47594
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Ankyrin B
Concentration:	1 mg/ml
Gene Name:	ANK2
Protein Name:	Ankyrin B
Other Name:	Ankyrin-2 (ANK-2) (Ankyrin-B) (Brain ankyrin) (Non-erythroid ankyrin)
Human Swiss Prot No.:	Q01484
Formulation:	Liquid in PBS containing 50% glycerol, and 0.24% sodium azide.
Specificity:	This antibody detects endogenous levels of Ankyrin B at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year