

Antibodies Product Datasheet

Product Name:	FSTL1Antibody
Catalog Number:	C47598
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FSTL1
Concentration:	1 mg/ml
Gene Name:	FSTL1 FRP
Protein Name:	FSTL1
Other Name:	Follistatin-related protein 1 (Follistatin-like protein 1)
Human Swiss Prot No.:	Q12841
Formulation:	Liquid in PBS containing 50% glycerol, and 0.29% sodium azide.
Specificity:	This antibody detects endogenous levels of FSTL1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year