

Antibodies Product Datasheet

Product Name:	NDR1/2Antibody
Catalog Number:	C47606
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDR1/2
Concentration:	1 mg/ml
Gene Name:	STK38 NDR1
Protein Name:	NDR1/2
Other Name:	Serine/threonine-protein kinase 38 (EC 2.7.11.1) (NDR1 protein kinase) (Nuclear Dbf2-related kinase 1)
Human Swiss Prot No.:	Q15208
Formulation:	Liquid in PBS containing 50% glycerol, and 0.37% sodium azide.
Specificity:	This antibody detects endogenous levels of NDR1/2 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year