

Antibodies Product Datasheet

Product Name:	GRK1/2Antibody
Catalog Number:	C47620
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GRK1/2
Concentration:	1 mg/ml
Gene Name:	ADRBK1 BARK BARK1 GRK2
Protein Name:	GRK1/2
Other Name:	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein coupled receptor kinase 2)
Human Swiss Prot No.:	P25098
Formulation:	Liquid in PBS containing 50% glycerol, and 0.51% sodium azide.
Specificity:	This antibody detects endogenous levels of GRK1/2 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year