

Antibodies Product Datasheet

Product Name:	HIVEP2 Antibody
Catalog Number:	C47623
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HIVEP2
Concentration:	1 mg/ml
Gene Name:	HIVEP2
Protein Name:	HIVEP2
Other Name:	Transcription factor HIVEP2 (Human immunodeficiency virus type I enhancer-binding protein 2) (HIV-EP2) (MHC-binding protein 2) (MBP-2)
Human Swiss Prot No.:	P31629
Formulation:	Liquid in PBS containing 50% glycerol, and 0.54% sodium azide.
Specificity:	This antibody detects endogenous levels of HIVEP2 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year