

Antibodies Product Datasheet

Product Name:	ISK7 Antibody
Catalog Number:	C47624
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ISK7
Concentration:	1 mg/ml
Gene Name:	SPINK7 ECG2 UNQ745/PRO1474
Protein Name:	ISK7
Other Name:	Serine protease inhibitor Kazal-type 7 (Esophagus cancer-related gene 2 protein) (ECRG-2)
Human Swiss Prot No.:	P58062
Formulation:	Liquid in PBS containing 50% glycerol, and 0.55% sodium azide.
Specificity:	This antibody detects endogenous levels of ISK7 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year