

Antibodies Product Datasheet

Product Name:	KITH_HHV11 Antibody
Catalog Number:	C47625
Reactivity:	HHV11
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KITH_HHV11
Concentration:	1 mg/ml
Gene Name:	#N/A
Protein Name:	KITH_HHV11
Other Name:	
Human Swiss Prot No.:	P0DTH5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.56% sodium azide.
Specificity:	This antibody detects endogenous levels of KITH_HHV11 at HHV11
Storage and Stability:	-20°C/1 year