

Antibodies Product Datasheet

Product Name:	HIF1AN/FIH-1 Antibody
Catalog Number:	C47634
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HIF1AN/FIH-1
Concentration:	1 mg/ml
Gene Name:	HIF1AN FIH1
Protein Name:	HIF1AN/FIH-1
Other Name:	Hypoxia-inducible factor 1-alpha inhibitor (EC 1.14.11.30) (EC 1.14.11.n4) (Factor inhibiting HIF-1) (FIH-1) (Hypoxia-inducible factor asparagine hydroxylase)
Human Swiss Prot No.:	Q9NWT6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.65% sodium azide.
Specificity:	This antibody detects endogenous levels of HIF1AN/FIH-1 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year