

Antibodies Product Datasheet

Product Name:	ITPR3 Antibody
Catalog Number:	C47638
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITPR3
Concentration:	1 mg/ml
Gene Name:	ITPR3
Protein Name:	ITPR3
Other Name:	Inositol 1,4,5-trisphosphate receptor type 3 (IP3 receptor isoform 3) (IP3R 3) (InsP3R3) (Type 3 inositol 1,4,5-trisphosphate receptor) (Type 3 InsP3 receptor)
Human Swiss Prot No.:	Q14573
Formulation:	Liquid in PBS containing 50% glycerol, and 0.69% sodium azide.
Specificity:	This antibody detects endogenous levels of ITPR3 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year