

Antibodies Product Datasheet

Product Name:	FCER1A Antibody
Catalog Number:	C47639
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FCER1A
Concentration:	1 mg/ml
Gene Name:	FCER1A FCE1A
Protein Name:	FCER1A
Other Name:	High affinity immunoglobulin epsilon receptor subunit alpha (Fc-epsilon RI-alpha) (FcERI) (IgE Fc receptor subunit alpha)
Human Swiss Prot No.:	P12319
Formulation:	Liquid in PBS containing 50% glycerol, and 0.70% sodium azide.
Specificity:	This antibody detects endogenous levels of FCER1A at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year