

Antibodies Product Datasheet

Product Name:	NMDARA2 Antibody
Catalog Number:	C47640
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NMDARA2
Concentration:	1 mg/ml
Gene Name:	GRIN2A NMDAR2A
Protein Name:	NMDARA2
Other Name:	Glutamate [NMDA] receptor subunit epsilon-1 (N-methyl D-aspartate receptor subtype 2A) (NMDAR2A) (NR2A) (hNR2A)
Human Swiss Prot No.:	Q12879
Formulation:	Liquid in PBS containing 50% glycerol, and 0.71% sodium azide.
Specificity:	This antibody detects endogenous levels of NMDARA2 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year