

Antibodies Product Datasheet

Product Name:	NDUFAF6 Antibody
Catalog Number:	C47642
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDUFAF6
Concentration:	1 mg/ml
Gene Name:	NDUFAF6 C8orf38
Protein Name:	NDUFAF6
Other Name:	NADH dehydrogenase (ubiquinone) complex I, assembly factor 6 (Putative phytoene synthase)
Human Swiss Prot No.:	Q330K2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.73% sodium azide.
Specificity:	This antibody detects endogenous levels of NDUFAF6 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year