

Antibodies Product Datasheet

Product Name:	glrx5 Antibody
Catalog Number:	C47645
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human glrx5
Concentration:	1 mg/ml
Gene Name:	GLRX5 C14orf87
Protein Name:	glrx5
Other Name:	Glutaredoxin-related protein 5, mitochondrial (Monothiol glutaredoxin-5)
Human Swiss Prot No.:	Q86SX6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.76% sodium azide.
Specificity:	This antibody detects endogenous levels of glrx5 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year