

Antibodies Product Datasheet

Product Name:	fxr1 Antibody
Catalog Number:	C47648
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human fxr1
Concentration:	1 mg/ml
Gene Name:	FXR1
Protein Name:	fxr1
Other Name:	Fragile X mental retardation syndrome-related protein 1 (hFXR1p)
Human Swiss Prot No.:	P51114
Formulation:	Liquid in PBS containing 50% glycerol, and 0.79% sodium azide.
Specificity:	This antibody detects endogenous levels of fxr1 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year