

Antibodies Product Datasheet

Product Name:	ACSL4 Antibody
Catalog Number:	C47650
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ACSL4
Concentration:	1 mg/ml
Gene Name:	ACSL4 ACS4 FACL4 LACS4
Protein Name:	ACSL4
Other Name:	Long-chain-fatty-acid--CoA ligase 4 (EC 6.2.1.3) (Long-chain acyl-CoA synthetase 4) (LACS 4)
Human Swiss Prot No.:	O60488
Formulation:	Liquid in PBS containing 50% glycerol, and 0.81% sodium azide.
Specificity:	This antibody detects endogenous levels of ACSL4 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year