

## Antibodies Product Datasheet

|                               |                                                                                           |
|-------------------------------|-------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | ApoL1 Monoclonal Antibody                                                                 |
| <b>Catalog Number:</b>        | C47662                                                                                    |
| <b>Reactivity:</b>            | Human                                                                                     |
| <b>Applications:</b>          | WB, ELISA                                                                                 |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.       |
| <b>Modification:</b>          | Unmodified/Total                                                                          |
| <b>Source:</b>                | Polyclonal Mouse                                                                          |
| <b>Purification:</b>          | Affinity purification                                                                     |
| <b>Immunogen:</b>             | Purified recombinant fragment of ApoL1 expressed in E. Coli.                              |
| <b>Concentration:</b>         |                                                                                           |
| <b>Gene Name:</b>             | APOL1                                                                                     |
| <b>Protein Name:</b>          | Apolipoprotein L1                                                                         |
| <b>Other Name:</b>            | APOL1; APOL; Apolipoprotein L1; Apolipoprotein L; Apo-L; ApoL; Apolipoprotein L-I; ApoL-I |
| <b>Human Swiss Prot No.:</b>  | O14791                                                                                    |
| <b>Formulation:</b>           | Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.                      |
| <b>Specificity:</b>           | ApoL1 Monoclonal Antibody detects endogenous levels of ApoL1 protein.                     |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                              |