

Antibodies Product Datasheet

Product Name:	IRAK4(p-T345/S346) Antibody
Catalog Number:	C47664
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IRAK4(p-T345/S346)
Concentration:	1 mg/ml
Gene Name:	IRAK4
Protein Name:	IRAK4(p-T345/S346)
Other Name:	Interleukin-1 receptor-associated kinase 4 (IRAK-4) (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-64)
Human Swiss Prot No.:	Q9NWZ3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.140% sodium azide.
Specificity:	This antibody detects endogenous levels of IRAK4(p-T345/S346) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year