

Antibodies Product Datasheet

Product Name:	HP1 gamma/CBX3 (citrulline R108) Antibody
Catalog Number:	C47668
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HP1 gamma/CBX3 (citrulline R108)
Concentration:	1 mg/ml
Gene Name:	CBX3
Protein Name:	HP1 gamma/CBX3 (citrulline R108)
Other Name:	Chromobox protein homolog 3 (HECH) (Heterochromatin protein 1 homolog gamma) (HP1 gamma) (Modifier 2 protein)
Human Swiss Prot No.:	Q13185
Formulation:	Liquid in PBS containing 50% glycerol, and 0.230% sodium azide.
Specificity:	This antibody detects endogenous levels of HP1 gamma/CBX3 (citrulline R108) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year