

Antibodies Product Datasheet

Product Name:	IFT88 Antibody
Catalog Number:	C48001
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IFT88
Concentration:	1 mg/ml
Gene Name:	IFT88 TG737 TTC10
Protein Name:	Intraflagellar transport protein 88 homolog (Recessive polycystic kidney disease protein Tg737 homolog) (Tetratricopeptide repeat protein 10) (TPR repeat protein 10)
Other Name:	
Human Swiss Prot No.:	Q13099
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IFT88 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year