

Antibodies Product Datasheet

Product Name:	FUS/TLS Antibody
Catalog Number:	C48003
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FUS/TLS
Concentration:	1 mg/ml
Gene Name:	FUS TLS
Protein Name:	RNA-binding protein FUS (75 kDa DNA-pairing protein) (Oncogene FUS) (Oncogene TLS) (POMp75) (Translocated in liposarcoma protein)
Other Name:	
Human Swiss Prot No.:	P35637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FUS/TLS at Human, Mouse
Storage and Stability:	Store at -20° C/1 year