

Antibodies Product Datasheet

Product Name:	CNOT1 Antibody
Catalog Number:	C48011
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CNOT1
Concentration:	1 mg/ml
Gene Name:	CNOT1 CDC39 KIAA1007 NOT1 AD-005
Protein Name:	CCR4-NOT transcription complex subunit 1 (CCR4-associated factor 1) (Negative regulator of transcription subunit 1 homolog) (NOT1H) (hNOT1)
Other Name:	
Human Swiss Prot No.:	A5YKK6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CNOT1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year