

Antibodies Product Datasheet

Product Name:	CAPRIN1 Antibody
Catalog Number:	C48018
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CAPRIN1
Concentration:	1 mg/ml
Gene Name:	CAPRIN1 GPIAP1 GPIP137 M11S1 RNG105
Protein Name:	Caprin-1 (Cell cycle-associated protein 1) (Cytoplasmic activation- and proliferation-associated protein 1) (GPI-anchored membrane protein 1) (GPI-anchored protein p137) (GPI-p137) (p137GPI) (Membrane component chromosome 11 surface marker 1) (RNA granule
Other Name:	
Human Swiss Prot No.:	Q14444
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CAPRIN1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year