

Antibodies Product Datasheet

Product Name:	IMPDH2 Antibody
Catalog Number:	C48027
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IMPDH2
Concentration:	1 mg/ml
Gene Name:	IMPDH2 IMPD2
Protein Name:	Inosine-5'-monophosphate dehydrogenase 2 (IMP dehydrogenase 2) (IMPD 2) (IMPDH 2) (EC 1.1.1.205) (IMPDH-II)
Other Name:	
Human Swiss Prot No.:	P12268
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IMPDH2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year