

Antibodies Product Datasheet

Product Name:	CHMP2A Antibody
Catalog Number:	C48030
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CHMP2A
Concentration:	1 mg/ml
Gene Name:	CHMP2A BC2 CHMP2
Protein Name:	Charged multivesicular body protein 2a (Chromatin-modifying protein 2a) (CHMP2a) (Putative breast adenocarcinoma marker BC-2) (Vacuolar protein sorting-associated protein 2-1) (Vps2-1) (hVps2-1)
Other Name:	
Human Swiss Prot No.:	O43633
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHMP2A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year