

Antibodies Product Datasheet

Product Name:	HRD1/SYVN1 Antibody
Catalog Number:	C48040
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HRD1/SYVN1
Concentration:	1 mg/ml
Gene Name:	SYVN1 HRD1 KIAA1810
Protein Name:	E3 ubiquitin-protein ligase synoviolin (EC 6.3.2.-) (Synovial apoptosis inhibitor 1)
Other Name:	
Human Swiss Prot No.:	Q86TM6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HRD1/SYVN1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year