

Antibodies Product Datasheet

Product Name:	CHCHD2 Antibody
Catalog Number:	C48043
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CHCHD2
Concentration:	1 mg/ml
Gene Name:	CHCHD2 C7orf17 AAG10
Protein Name:	Coiled-coil-helix-coiled-coil-helix domain-containing protein 2, mitochondrial (Aging-associated gene 10 protein) (HCV NS2 trans-regulated protein) (NS2TP)
Other Name:	
Human Swiss Prot No.:	Q9Y6H1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHCHD2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year