

Antibodies Product Datasheet

Product Name:	DHX9 Antibody
Catalog Number:	C48060
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHX9
Concentration:	1 mg/ml
Gene Name:	DHX9 DDX9 LKP NDH2
Protein Name:	ATP-dependent RNA helicase A (RHA) (EC 3.6.4.13) (DEAH box protein 9) (Leukophysin) (LKP) (Nuclear DNA helicase II) (NDH II)
Other Name:	
Human Swiss Prot No.:	Q08211
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHX9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year