

Antibodies Product Datasheet

Product Name:	IFT57 Antibody
Catalog Number:	C48076
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IFT57
Concentration:	1 mg/ml
Gene Name:	IFT57 DERP8 ESRRBL1 HIPPI
Protein Name:	Intraflagellar transport protein 57 homolog (Dermal papilla-derived protein 8) (Estrogen-related receptor beta-like protein 1) (HIP1-interacting protein) (MHS4R2)
Other Name:	
Human Swiss Prot No.:	Q9NWB7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IFT57 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year