

Antibodies Product Datasheet

Product Name:	AHCYL1 Antibody
Catalog Number:	C48077
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AHCYL1
Concentration:	1 mg/ml
Gene Name:	AHCYL1 DCAL XPVKONA
Protein Name:	Putative adenosylhomocysteinase 2 (AdoHcyase 2) (EC 3.3.1.1) (DC-expressed AHCY-like molecule) (S-adenosyl-L-homocysteine hydrolase 2) (S-adenosylhomocysteine hydrolase-like protein 1)
Other Name:	
Human Swiss Prot No.:	O43865
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AHCYL1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year