

Antibodies Product Datasheet

Product Name:	ASF/SF2 Antibody
Catalog Number:	C48079
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASF/SF2
Concentration:	1 mg/ml
Gene Name:	SRSF1 ASF SF2 SF2P33 SFRS1 OK/SW-cl.3
Protein Name:	Serine/arginine-rich splicing factor 1 (Alternative-splicing factor 1) (ASF-1) (Splicing factor, arginine/serine-rich 1) (pre-mRNA-splicing factor SF2, P33 subunit)
Other Name:	
Human Swiss Prot No.:	Q07955
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASF/SF2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year