

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                     |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Kindlin 2 Antibody                                                                                                                                                                  |
| <b>Catalog Number:</b>        | C48080                                                                                                                                                                              |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                                                                                        |
| <b>Applications:</b>          | WB                                                                                                                                                                                  |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                                                                                       |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                    |
| <b>Source:</b>                | Polyclonal, Rabbit, IgG                                                                                                                                                             |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                               |
| <b>Immunogen:</b>             | Synthesized peptide derived from human Kindlin 2                                                                                                                                    |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                             |
| <b>Gene Name:</b>             | FERMT2 KIND2 MIG2 PLEKHC1                                                                                                                                                           |
| <b>Protein Name:</b>          | Fermitin family homolog 2 (Kindlin-2) (Mitogen-inducible gene 2 protein) (MIG-2) (Pleckstrin homology domain-containing family C member 1) (PH domain-containing family C member 1) |
| <b>Other Name:</b>            |                                                                                                                                                                                     |
| <b>Human Swiss Prot No.:</b>  | Q96AC1                                                                                                                                                                              |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                             |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Kindlin 2 at Human, Mouse                                                                                                                |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                                                                              |