

Antibodies Product Datasheet

Product Name:	JTV1 Antibody
Catalog Number:	C48083
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human JTV1
Concentration:	1 mg/ml
Gene Name:	AIMP2 JTV1 PRO0992
Protein Name:	Aminoacyl tRNA synthase complex-interacting multifunctional protein 2 (Multisynthase complex auxiliary component p38) (Protein JTV-1)
Other Name:	
Human Swiss Prot No.:	Q13155
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of JTV1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year