

Antibodies Product Datasheet

Product Name:	INO80 Antibody
Catalog Number:	C48102
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human INO80
Concentration:	1 mg/ml
Gene Name:	INO80 INO80A INOC1 KIAA1259
Protein Name:	DNA helicase INO80 (hINO80) (EC 3.6.4.12) (INO80 complex subunit A) (Putative DNA helicase INO80 complex homolog 1)
Other Name:	
Human Swiss Prot No.:	Q9ULG1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of INO80 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year