

Antibodies Product Datasheet

Product Name:	FATP2 Antibody
Catalog Number:	C48108
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FATP2
Concentration:	1 mg/ml
Gene Name:	SLC27A2 ACSVL1 FACVL1 FATP2 VLACS
Protein Name:	Very long-chain acyl-CoA synthetase (VLACS) (VLCS) (EC 6.2.1.-) (Fatty acid transport protein 2) (FATP-2) (Fatty-acid-coenzyme A ligase, very long-chain 1) (Long-chain-fatty-acid--CoA ligase) (EC 6.2.1.3) (Solute carrier family 27 member 2) (THCA-CoA liga
Other Name:	
Human Swiss Prot No.:	O14975
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FATP2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year