

Antibodies Product Datasheet

Product Name:	IFT25 Antibody
Catalog Number:	C48119
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IFT25
Concentration:	1 mg/ml
Gene Name:	HSPB11 C1orf41 HSPC034
Protein Name:	Heat shock protein beta-11 (Hspb11) (Placental protein 25) (PP25)
Other Name:	
Human Swiss Prot No.:	Q9Y547
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IFT25 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year