

Antibodies Product Datasheet

Product Name:	DHRS3 Antibody
Catalog Number:	C48126
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHRS3
Concentration:	1 mg/ml
Gene Name:	DHRS3 UNQ2424/PRO4983
Protein Name:	Short-chain dehydrogenase/reductase 3 (EC 1.1.1.300) (DD83.1) (Retinal short-chain dehydrogenase/reductase 1) (retSDR1)
Other Name:	
Human Swiss Prot No.:	O75911
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHRS3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year